



XPERT® NOROVIRUS

Test Reagent Kit	Xpert® Norovirus
Technology	Real time PCR
Targets (unique to Cepheid)	GI, GII
Batch or On Demand	On Demand
Minimum Batch Size	1
Sample Extraction	Automated/Integrated
Precision Pipetting	Not Required
TAT	Results in as little as one hour*
Hands-On Time	<1 minute
Limit of Detection (LoD) For Common Strains	GI 1.87 x 10e5 copies/mL GII 1.55 x 10e4 copies/mL
Internal Controls: Process (SPC)	✓
Internal Controls: Probe Function/Detection (PCC)	✓
Sensitivity	GI 100% GII 100%
Specificity	GI 99.5% GII 98.9%
Sample Stability	2 days at 2-8°C
Kit Storage	2-8°C
Commercial Controls	Contact Cepheid Technical Support

* with EAT (Early Assay Termination), for positive results



In Vitro Diagnostic Medical Device

CORPORATE HEADQUARTERS

904 Caribbean Drive
Sunnyvale, CA 94089 USA

TOLL FREE 1.888.336.2743
PHONE 1.408.541.4191
FAX 1.408.541.4192

EUROPEAN HEADQUARTERS

Vira Solelh
81470 Maurens-Scopont France

PHONE 33.563.82.53.00
FAX 33.563.82.53.01
EMAIL cepheid@cepheideurope.fr

www.Cepheidinternational.com